

TEST KIT COMPARISON

FUELSTAT®*Plus*



FUELSTAT®*One*



Immunoassay antibody
detection rapid lateral flow test

Six competitive
lateral flow devices

Detects **specific** fuel degrading
micro-organisms

Semi-quantitative result:
Negligible / **Moderate** / **Heavy**

Result in **15-20 minutes**

Result can be read by eye,
verified using free FUELSTAT®
Result app, instant PDF report
and FUELSTAT® *Report* portal

Available as single test kits
or in boxes of 8

Protein-ligand complex
detection rapid lateral flow test

Single sandwich
lateral flow device

Detects a **broad range** of
micro-organisms associated
with fuel tanks

Fully-quantitative result:
Negligible / **Moderate** / **Heavy**
with position indicator and
numerical contamination value

Result in **20-30 minutes**

Result read using free
FUELSTAT® *Result* app, instant
PDF report and FUELSTAT®
Report portal

Available as single test kits
or in boxes of 8